# Overview Présentation / Überblick / Panoramica / Overzicht / Opis ogólny / Обзор / Vista general

1

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
Figure 3 displays the typical spectral transmittance curves for YP2.2 / CP2.2 with a dioptric power of +20.0 together with the spectral transmittance curve\* for the phakic eye of a 53-year-old patient.

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












**PRESCRIPTION DEVICE**

Caution: This device is restricted to sale by or on the order of an ophthalmologist.

**DEVICE DESCRIPTION**

The Avansee™ Preload1P / Avansee™ Preload1P Clear from *Kowa Company, Ltd.* is a foldable posterior chamber aspheric intraocular lens (IOL) (Fig.2) that is preloaded in a single-use *Kowa* original injector (Fig.1). The optic of the lens and the **modified-C** haptics are made from a UV-absorbing hydrophobic soft acrylic material which, in the case of the Yellow type also contains proprietary blue-light filtering. The spectral transmittance of the IOL (Yellow type) closely replicates that of the natural crystalline lens (Fig.3). Prior to insertion, the optical portion of the lens is folded to allow for placement through a small incision. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance.

**INTENDED USE**

*Kowa's* Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This allows the lens to function as a refractive medium in the visual correction of aphakia.

**INDICATIONS**

*Kowa's* Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is placed in a capsular bag and is designed for implantation after extracapsular cataract extraction or phacoemulsification of cataracts.

**WARNINGS**

1. Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting the Avansee™ Preload1P / Avansee™ Preload1P Clear lens in a child or a "special care patient" who has one or more of the following conditions. Before surgery Corneal endothelial damage / Glaucoma / Uveitis / Diabetic retinopathy / Retinal detachment / Congenital ocular anomalies / Choroidal haemorrhage / Shallow anterior chamber / Microphthalmos / Corneal dystrophy / Optic nerve atrophy / Ocular hypertension / Amydriasis / Amblyopia / History of keratoplasty / Iritis / Corneal anomalies / Macular degeneration / Retinal degeneration / Atopic disease / Pseudoexfoliation syndrome and zonular weakness / Zonular rupture and lens luxation (including lens subluxation) / Rubeosis iridis / Patients with intraoperative occurrences of any severe adverse event / or patients who, assessed by a surgeon, require special care for reasons such as an accompanying systemic or ocular disease.
2. Prior to surgery, the risks and benefits associated with the implantation of the IOL should be clearly explained by the surgeon to the patient.
3. Additional attention is required for "special care patients", including a post-operative follow up by an experienced ophthalmologist with adequate equipment, due to the higher risk of complications and/or insufficient vision recovery. Implantation of the IOL in children should be performed by an ophthalmologist who has sufficient knowledge and experience in paediatric care. This is particularly important when treating a child (under 2 years old) whose small eyeballs can make implantation and the handling of instruments difficult. The chances of IOL exchange are also higher due to the change of eye axial length with the child's growth. Therefore it is important to provide the parents of the young patient with informed consent information prior to any surgery and the IOL should be implanted to the young patients carefully after considering the benefit/risk ratio.
4. For active uveitis and a child with uveitis, inflammation should be suppressed by medical treatment prior to intraocular lens implantation, as surgical invasion may cause aggravation of the uveitis or another complication.
5. As with any surgical procedure, potential complications accompanying IOL implantation can occur. Adverse events and malfunctions may require discontinuation of the implantation, since they can lead to blindness, permanent problems with visual acuity and grave health hazard or may require IOL extraction and IOL exchange.
6. Complications accompanying IOL implantation may include but are not limited to:  
<Adverse events>  
Corneal oedema / Keratitis (including corneal erosion) / Corneal endothelial damage / Acute corneal decompensation / Detachment of Descemet's membrane / Conjunctivitis and subconjunctival haemorrhage / Hyphaema / Hypopyon / Iris damage / Iritis (iridocyclitis) / Iris adhesion / Iris prolapse / Abnormal pupil (block, capture, deformation, dilatation, etc.) / Uveitis / Zonular rupture / Posterior capsular rupture / After cataract / Vitreous haemorrhage and opacity / Vitreous prolapse / Detachment, hole, tear, etc. of retinal tissue (including macula) / Retinal detachment / Choroidal detachment / Choroid haemorrhage / Macular oedema and degeneration / Expulsive haemorrhage / Endophthalmitis / Fibrin precipitation /

Secondary glaucoma / Intraocular pressure elevation (including transient elevated intraocular pressure and ocular hypertension) / Intraocular pressure lowering / Dysphotopsia / Deterioration of visual function (visual acuity and contrast sensitivity) / Error of predicted refractive power / Wound leak.  
<Malfunctions>

- <Malfunctions>  
Optic damage (breakage, scratch, etc.) / Haptic damage (breakage, scratch, detachment, etc.) / Adhesion of foreign bodies on the lens surface / Lens surface reflection / Lens discoloration or pseudocoloration / Lens opacification (including glistening) / Lens luxation / Lens decentration / Lens dislocation / Lens jamming.
7. Back-ups for the intraocular lens and the insertion device should be made available prior to surgery as a part of emergency preparations.
  8. DO NOT re-use the lens. For single use only.
  9. DO NOT re-sterilise. Re-use or re-sterilisation of any component may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Re-use or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to the transmission of infectious disease(s) from one patient to another.
  10. DO NOT store the lens in direct sunlight or in any hot and humid places. Store the lens at room temperature (between 1 °C and 30 °C) and keep it dry.

**PRECAUTIONS**

1. Before using this device, please ensure that you have carefully read and understood this document in order to complete the procedure safely.
2. This device should be placed in the posterior chamber. Implantation of the lens in the anterior chamber has not been evaluated for safety and efficacy.
3. Before opening the Avansee™ Preload1P / Avansee™ Preload1P Clear packaging, confirm the description including the IOL model, lens power, and expiration date located on the label.
4. The temperatures of the device and the ophthalmic viscosurgical device (OVD) during use should be kept between 21 °C and 25 °C.
5. The inside of the aluminium packaging is not sterile. In a non-sterile environment, remove its contents including the primary (blister) package from the aluminium packaging.
6. Should any damage or abnormality be found on the primary package after opening the aluminium packaging, DO NOT use, as the device may no longer be sterile.
7. After opening the primary package, the device should be handled aseptically.
8. When removing the device from the primary package, ensure to be careful not to drop it.
9. Before using the device, carefully examine all parts for damage, or any other abnormality. Should any abnormality be found, DO NOT use it.
10. If the lens stage is dislodged before use, DO NOT use it.
11. DO NOT remove the lens stage before the injection of the OVD.
12. DO NOT open the cap on the injector body.
13. If there are any abnormalities during the procedure, immediately stop using this device.
14. When implantation has been stopped during the procedure due to an abnormality, DO NOT re-use the device and discard it.

**INSTRUCTIONS FOR USE**

<Preparation for lens implantation>

1. Prior to implanting, examine the lens package for type, power and proper configuration.
2. Open the primary (blister) package and remove the device in a sterile environment.
3. Insert the injection needle for the OVD into the inlet located in the cap on the injector body, and fill the nozzle close to the dashed line shown in the diagram (Fig. 4) with the OVD in order to inject at least approx. 0.17 ml.

<Lens implantation>

4. With the *Kowa* mark portion of the injector body firmly in hand, hold the lens stage between your thumb and middle finger on both sides and support its fore part with your index finger (Fig. 5), and then remove the lens stage slowly away from the injector body towards the nozzle tip (Fig. 6).
5. Push the plunger at a slow and constant rate; move the IOL forward and stop it at the point where the IOL optic is rolled and its edges make secure contact (Fig. 7). Once this has been done, immediately (within 20 seconds) place the IOL into the eye.
6. Insert the nozzle tip with the bevel down through an incision to just before the central pupillary area (Fig. 8).
7. While pushing the plunger ahead at a constant rate release the IOL inside the capsular bag.

- Push until the trailing haptic is completely extruded and the IOL is released, and pull the injector body out of the anterior chamber. Then, stabilise the IOL inside the capsular bag using a hook or other tools.
- Care should be taken to remove the OVD from the eye at the end of surgery.
- Cautions in handling:**
  - A sufficient amount of OVD should be injected to fill the nozzle and cover the entire lens optic by inserting the injection needle into the inlet located in the cap on the injector body (to the level of at least 0.17 ml from the inlet). Insufficient OVD may damage the IOL because of inadequate lubrication inside the nozzle.
  - DO NOT leave the injector with the OVD filled inside the nozzle. Immediately start the next procedure.
  - DO NOT advance the plunger before removing the lens stage.
  - When removing the lens stage, DO NOT twist it laterally (Fig. 9).
  - Upon removing the lens stage, confirm that there is no abnormality in the IOL. Should a haptic of the IOL become deformed and protrude out of the hole in the back of the injector body as shown in Fig. 10, DO NOT use it.
  - After removing the lens stage, DO NOT fill it with additional OVD.
  - DO NOT push the plunger rapidly or intermittently. Push the plunger forward at a constant rate.
  - DO NOT pull back the plunger once the plunger has been pushed forward.
  - Keep the bevel down while pushing the plunger ahead at a constant rate to release the IOL inside the capsular bag.
  - If you feel strong resistance or any abnormal feeling during advancement of the plunger, stop using this device.
  - During advancement of the plunger, stop using this device if a trailing haptic is out of the groove along the plunger, or if :
    - the leading haptic should twist and extend forward (Fig. 11);
    - the leading haptic should bend and stretch out (Fig. 12);
    - the trailing haptic should extend out (Fig. 13);
    - the plunger should pass above or under the lens optic or bend the optic irregularly (Fig. 14); or
    - the plunger should go too far toward the right or left side (Fig. 15).
  - The IOL should be moved forward until the lens optic is rolled and its edges make secure contact. If you stop advancement short of this location, the plunger may get into and be caught in the lens optic and the lens may be damaged when pushing and inserting the IOL afterward.

- After the IOL is moved into the near tip of the nozzle, move on to the next procedure promptly (within 20 seconds) and place the lens into the eye. Failure of prompt implantation may cause the lens to become adhered inside the nozzle, possibly resulting in damage to the IOL or improper implantation.
- During the IOL implantation procedure, DO NOT push the lens to the extent of being exposed from the nozzle tip (Fig. 16).
- The IOL should be released slowly after the nozzle tip has been located just before the centre of the pupillary area.

**CALCULATION OF LENS POWER**

The surgeon should determine preoperatively the power of lens to be implanted. The A-constant below is presented as a reference value for lens power calculations. When calculating the exact lens power, it is recommended that calculations should be performed individually based on the equipment used and the operating surgeon's experience.

Model	A-constant (A-mode)
YP2.2 / CP2.2	118.6

Lens power calculation methods are described in the following references.

Sanders, D.R., Retzlaff, J., and Kraff, M.C., "Comparison of SRK II formula and other second generation formulae", *Journal of Cataract and Refractive Surgery*, Vol. 14, pp.136-141, 1988.  
Retzlaff, J.A., Sanders, D.R., and Kraff, M.C., "Development of the SRK/T intraocular lens implant power calculation formula", *Journal of Cataract and Refractive Surgery*, Vol. 16, pp.333-340, 1990; ERRATA, Vol. 16, pp.528, 1990.

**PATIENT IDENTIFICATION CARD**

The Patient Identification Card included in the outer packaging should be completed and given to the patient. The patient must provide this card when seeking any future eye care.

**HOW SUPPLIED**

The Avansee™ Preload1P / Avansee™ Preload1P Clear is supplied sterile in a primary (blister) package for single use only. The primary (blister) package is sterilised with ethylene oxide gas and the package should be opened only under sterile conditions. DO NOT re-sterilise.

**EXPIRATION DATE (sterility expiration date)**

The sterility expiration date is clearly indicated on the outer package. The date is expressed as four digits for the year and two digits for the month.

**RETURN/EXCHANGE POLICY**

Regarding lens return or exchange, please contact with your local distributor or Kowa Pharmaceutical Europe Co. Ltd. located at 105 Wharfedale Road, Winnersh Triangle, Wokingham, Berkshire, RG41 5RB, U.K.

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Symbols	Description	Symbols	Description	Symbols	Description
	Do not re-use / For single use only		Caution		Manufacturer
	Do not re-sterilise		Keep away from sunlight		Authorized representative in the European Community
	Do not use if package is damaged		Keep dry		Dioptre
	Sterilised using ethylene oxide		Temperature limitation	$\phi_T$	Overall diameter
	Non-sterile		Serial number	$\phi_B$	Optic diameter (Body diameter)
	Consult instructions for use		Use by (YYYY-MM : year-month)		

Manufacturer:  
Kowa Company, Ltd.  
4-14, 3-Chome, Nihonbashi-Honcho, Chuo-ku, Tokyo 103-8433, JAPAN  
Authorized Representative in the European Community:  
Kowa Pharmaceutical Europe Co. Ltd.  
105 Wharfedale Road, Winnersh Triangle, Wokingham, Berkshire, RG41 5RB, U.K.  
Phone: 00-800-2826-7333  
Fax: +44 (0) 118-944-3806  
www.kowapharmaceuticals.eu  
www.avansee.eu

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Kowa Pharmaceutical Europe Co. Ltd.

105 Wharfedale Road, Winkersham Triangle, Wokingham, Berkshire, RG41 5RB, U.K.

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Sanders, D.R., Retzlaff, J., and Kraff, M.C., "Comparison of SRK II formula and other second generation formulae", *Journal of Cataract and Refractive Surgery*, Vol. 14, pp.136-141, 1988.

Retzlaff, J.A., Sanders, D.R., and Kraff, M.C., "Development of the SRK/T intraocular lens implant power calculation formula", *Journal of Cataract and Refractive Surgery*, Vol. 16, pp.333-340, 1990; ERRATA, Vol. 16, pp.528, 1990.

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