

X Do Not


À ne pas faire / Verboden / Divieti / Nooit / Przeciwwskazania / NA MHN / No


- 9** DO NOT twist laterally when removing the lens stage.
NE PAS provoquer de torsion latérale lors du retrait du support de stockage de la lentille.
Linsenhalter beim Entfernen NICHT seitlich drehen.
NON ruotare lateralmente durante la rimozione della sicura della lente.
NIET zijaarst draaien bij verwijdering van de lenspatroon.
NIE odginać na bok podczas zdejmowania osłony soczewki.
NA MHN συστρέφετε πλευρικά τη θήκη του φακού κατά την αφαίρεση της.
NO torcer lateralmente al retirar la plataforma de la lente.


- 10** DO NOT use – should haptic become deformed or protrude from the device.
NE PAS utiliser – si l'haptique est déformé ou s'il dépasse de l'appareil.
NICHT verwenden – Haptik verformt oder ragt aus dem Produkt heraus.
NON utilizzare in caso di deformazione o di protrusione di un'apica dal dispositivo.
NIET gebruiken – als het haptische deel vervormd is of uit het instrument steekt.
NIE używać jeżeli część haptyczna uległa odkształceniu lub wystaje z urządzenia.
ΜΗΝ χρησιμοποιείτε – όταν το απτικό στοιχείο έχει παραμορφωθεί ή προεξέχει από τη συσκευή.
NO utilizar – si el háptico se deforma o sobresale del dispositivo.


- 11** DO NOT use – leading haptic twisted and/or extended forward.
NE PAS utiliser – si l'haptique antérieur est tordu ou s'étend vers l'avant.
NICHT verwenden – vordere Haptik verdreht und/oder nach vorne gestreckt.
NON utilizzare in caso di apica di testa ruotata e/o allungata in avanti.
NIET gebruiken – als het voorste haptische deel gedraaid is en/of naar voren steekt.
NIE używać – przednia część haptyczna skręcona i wyciągnięta ku przodowi.
ΜΗΝ χρησιμοποιείτε – όταν το άνω άκρο του απτικού στοιχείου έχει στρεβλωθεί και/ή εκτείνεται προς τα εμπρός.
NO utilizar – háptico delantero retorcido y/o extendido.


- 12** DO NOT use – leading haptic bent or stretched out.
NE PAS utiliser – si l'haptique antérieur est coulé ou dépasse.
NICHT verwenden – vordere Haptik verbogen oder nach außen gestreckt.
NON utilizzare – apica di testa piegata o distesa.
NIET gebruiken – voorste haptische deel verbogen of gerekt.
NIE używać: przednia część haptyczna zgięta lub rozciągnięta.
NA MHN χρησιμοποιείτε – όταν η πρόσθια ακύλη έχει λυγίσει ή εκτείνεται εκτός.
NO utilizar – háptico delantero doblado o extendido.


- 13** DO NOT use – trailing haptic extended out.
NE PAS utiliser – si l'haptique postérieur est détendu.
NICHT verwenden – hintere Haptik herausgestreckt.
NON utilizzare – apica di coda allungata verso l'esterno.
NIET gebruiken – achterste haptische deel steekt uit.
NIE używać: tylna część haptyczna rozciągnięta.
NA MHN χρησιμοποιείτε – όταν η οπίσθια ακύλη επεκτείνεται προς τα έξω.
NO utilizar – háptico trasero extendido.

- 14** DO NOT use – should the plunger pass above or under the lens optic or bend the optic irregularly.
NE PAS utiliser – si le piston entre dans le corps de la lentille ou s'il tord anormalement la lentille.
NICHT verwenden – Kolben über oder unter die Linsenoptik oder auf ungewöhnliche Weise in der Optik gebogen.
NON utilizzare – se lo stantuffo passa sopra o sotto l'ottica della lente o piega l'ottica in modo irregolare.
NIET gebruiken – indien de zuiger boven of onder de optische lens loopt of de optische lens onjuist buigt.
NIE używać, jeżeli tłok przechodzi powyżej lub poniżej części optycznej soczewki lub zagina tę część w nieregularny sposób.
NA MHN χρησιμοποιείτε – όταν το έμβολο έχει περάσει πάνω ή κάτω από την οπτική διάταξη του φακού ή έχει κάμψει ακανόνιστα το φακό.
NO utilizar – si el émbolo pasara por encima o por debajo de la óptica de la lente o doblara la óptica de manera irregular.

- 15** DO NOT use – plunger has moved too far towards the left or right side.
NE PAS utiliser – si le piston est allé trop à droite ou à gauche.
NICHT verwenden – Kolben zu weit nach links oder rechts verschoben.
NON utilizzare – stantuffo posizionato troppo a destra o a sinistra.
NIET gebruiken – zuiger te ver naar links of rechts bewegen.
NIE używać: tłok przesunął się zbyt daleko w lewo lub w prawo.
NA MHN χρησιμοποιείτε – όταν το έμβολο έχει προχωρήσει πολύ προς τη αριστερή ή τη δεξιά πλευρά.
NO utilizar – el émbolo se ha movido con demasiada rapidez hacia el lado izquierdo o derecho.

- 16** DO NOT use – lens exposed at nozzle tip before insertion.
NE PAS utiliser – si la lentille a été exposée à l'extrémité de l'embout avant insertion.
NICHT verwenden – Linse vor dem Einsetzen auf Injektorspitze gelangt.
NON utilizzare – lente sporgente sulla punta dell'ugello prima dell'inserimento.
NIET gebruiken – lens bij punt Injector opening blootgesteld voor het inbrengen.
NIE używać: soczewka wystaje z końcówki dyszy przed wprowadzeniem.
NA MHN χρησιμοποιείτε – όταν ο φακός έχει βγει από το άκρο του ρύγχους πριν την εισαγωγή στο μάτι.
NO utilizar – lente expuesta en la punta de la boquilla antes de la inserción.

PRESCRIPTION DEVICE

Caution: This device is restricted to sale by or on the order of an ophthalmologist.

DEVICE DESCRIPTION

The Avansee™ Preload1P / Avansee™ Preload1P Clear from Kowa Company, Ltd. is a foldable posterior chamber aspheric intraocular lens (IOL) (Fig.2) that is preloaded in a single-use Kowa original injector (Fig.1). The optic of the lens and the modified-C haptics are made from a UV-absorbing hydrophobic soft acrylic material which, in the case of the Yellow type also contains proprietary blue-light filtering. The spectral transmittance of the IOL (Yellow type) closely replicates that of the natural crystalline lens (Fig.3). Prior to insertion, the optical portion of the lens is folded to allow for placement through a small incision. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance.

INTENDED USE

Kowa's Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This allows the lens to function as a refractive medium in the visual correction of aphakia.

INDICATIONS

Kowa's Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is placed in a capsular bag and is designed for implantation after extracapsular cataract extraction or phacemulsification of cataracts.

WARNINGS

- Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting the Avansee™ Preload1P / Avansee™ Preload1P Clear lens in a child or a "special care patient" who has one or more of the following conditions. Before surgery Corneal endothelial damage / Glaucoma / Uveitis / Diabetic retinopathy / Retinal detachment / Congenital ocular anomalies / Choroidal haemorrhage / Shallow anterior chamber / Microphthalmos / Corneal dystrophy / Optic nerve atrophy / Ocular hypertension / Amyotriasis / Amblyopia / History of keratoplasty / Iritis / Corneal anomalies / Macular degeneration / Retinal degeneration / Atopic disease / Pseudoexfoliation syndrome and zonular weakness / Zonular rupture and lens luxation (including lens subluxation) / Rubeosis iridis / Patients with intraoperative occurrences of any severe adverse event / or patients who, assessed by a surgeon, require special care for reasons such as an accompanying systemic or ophthalmic disease.
- Prior to surgery, the risks and benefits associated with the implantation of the IOL should be clearly explained by the surgeon to the patient.
- Additional attention is required for "special care patients", including a post-operative follow up by an experienced ophthalmologist with adequate equipment, due to the higher risk of complications and/or insufficient vision recovery. Implantation of the IOL in children should be performed by an ophthalmologist who has sufficient knowledge and experience in paediatric care. This is particularly important when treating a child (under 2 years old) whose small eyeballs can make implantation and the handling of instruments difficult. The chances of IOL exchange are also higher due to the change of eye axial length with the child's growth. Therefore it is important to provide the parents of the young patient with informed consent information prior to any surgery and the IOL should be implanted to the young patients carefully after considering the benefit/risk ratio.
- For active uveitis and a child with uveitis, inflammation should be suppressed by medical treatment prior to intraocular lens implantation, as surgical invasion may cause aggravation of the uveitis or another complication.
- As with any surgical procedure, potential complications accompanying IOL implantation can occur. Adverse events and malfunctions may require discontinuation of the implantation, since they can lead to blindness, permanent problems with visual acuity and grave health hazard or may require IOL extraction and IOL exchange.
- Complications accompanying IOL implantation may include but are not limited to:
 - <Adverse events>
 - Corneal oedema / Keratitis (including corneal erosion) / Corneal endothelial damage / Acute corneal decompensation / Detachment of Descemet's membrane / Conjunctivitis and subconjunctival haemorrhage / Hyphaema / Hypopyon / Iris damage / Iritis (iridocyclitis) / Iris adhesion / Iris prolapse / Abnormal pupil (block, capture, deformation, dilatation, etc.) / Uveitis / Zonular rupture / Posterior capsular rupture / After cataract / Vitreous haemorrhage and opacity / Vitreous prolapse / Detachment, hole, tear, etc. of retinal tissue (including macula) / Retinal detachment / Choroidal detachment / Choroid haemorrhage / Macular oedema and degeneration / Expulsive haemorrhage / Endophthalmitis / Fibrin precipitation /

Secondary glaucoma / Intraocular pressure elevation (including transient elevated intraocular pressure and ocular hypertension) / Intraocular pressure lowering / Dysphotopsia / Deterioration of visual function (visual acuity and contrast sensitivity) / Error of predicted refractive power / Wound leak.

<Malfunctions>

- Optic damage (breakage, scratch, etc.) / Haptic damage (breakage, scratch, detachment, etc.) / Adhesion of foreign bodies on the lens surface / Lens surface reflection / Lens discoloration or pseudodiscoloration / Lens opacification (including glistening) / Lens luxation / Lens decentration / Lens dislocation / Lens jamming.
- Back-ups for the intraocular lens and the insertion device should be made available prior to surgery as a part of emergency preparations.
 - DO NOT re-use the lens. For single use only.
 - DO NOT re-sterilise. Re-use or re-sterilisation of any component may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Re-use or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to the transmission of infectious disease(s) from one patient to another.
 - DO NOT store the lens in direct sunlight or in any hot and humid places. Store the lens at 25 °C or lower and keep it dry.

PRECAUTIONS

- Before using this device, please ensure that you have carefully read and understood this document in order to complete the procedure safely.
- DO NOT use this lens for the patients with capsular rupture or with zonular rupture, or with large posterior capsulotomy.
- This device should be placed in the posterior chamber. Implantation of the lens in the anterior chamber has not been evaluated for safety and efficacy. The lens should not be placed in the ciliary sulcus.
- Before opening the Avansee™ Preload1P / Avansee™ Preload1P Clear packaging, confirm the description including the IOL model, lens power, and expiration date located on the label.
- The temperatures of the device and the ophthalmic viscosurgical device (OVD) during use should be kept between 21 °C and 25 °C.
- The inside of the aluminium packaging is not sterile. In a non-sterile environment, remove its contents including the primary (blister) package from the aluminium packaging.
- Should any damage or abnormality be found on the primary package after opening the aluminium packaging, DO NOT use, as the device may no longer be sterile.
- After opening the primary package, the device should be handled aseptically.
- When removing the device from the primary package, ensure to be careful not to drop it.
- Before using the device, carefully examine all parts for damage, or any other abnormality. Should any abnormality be found, DO NOT use it.
- If the lens stage is dislodged before use, DO NOT use it.
- DO NOT remove the lens stage before the injection of the OVD.
- DO NOT open the cap on the injector body.
- If there are any abnormalities during the procedure, immediately stop using this device.
- Except for aborting lens insertion, do not interrupt insertion after the lens starts coming out from the nozzle, and keep pushing the plunger until trailing haptic is completely extruded from the nozzle.
- When implantation has been stopped during the procedure due to an abnormality, DO NOT re-use the device and discard it.

INSTRUCTIONS FOR USE**<Preparation for lens implantation>**

- Prior to implanting, examine the lens package for type, power and proper configuration.
- Open the primary (blister) package and remove the device in a sterile environment.
- Insert the injection needle for the OVD into the inlet located in the cap on the injector body, and fill the nozzle close to the dashed line shown in the diagram (Fig. 4) with the OVD in order to inject at least approx. 0.17 ml.

<Lens implantation>

- With the Kowa mark portion of the injector body firmly in hand, hold the lens stage between your thumb and middle finger on both sides and support its fore part with your index finger (Fig. 5), and then remove the lens stage slowly away from the injector body straight (Fig. 6).
- Push the plunger at a slow and constant rate; move the IOL forward and stop it at the point where the IOL optic is rolled and its edges make secure contact (Fig. 7). Once this has been done, immediately (within 20 seconds) place

- the IOL into the eye.
- 6. Insert the nozzle tip with the inlet (*Kowa* mark) upward through an incision to just before the central pupillary area (Fig. 8).
- 7. While pushing the plunger ahead at a constant rate release the IOL inside the capsular bag.
- 8. Push until the trailing haptic is completely extruded and the IOL is released, and pull the injector body out of the anterior chamber. Then, stabilise the IOL inside the capsular bag using a hook or other tools.
- 9. Care should be taken to remove the OVD from the eye at the end of surgery.

Cautions in handling:

- 1. A sufficient amount of OVD should be injected to fill the nozzle and cover the entire lens optic by inserting the injection needle into the inlet located in the cap on the injector body (to the level of at least 0.17 ml from the inlet). Insufficient OVD may damage the IOL because of inadequate lubrication inside the nozzle.
- 2. DO NOT leave the injector with the OVD filled inside the nozzle. Immediately start the next procedure.
- 3. DO NOT advance the plunger before removing the lens stage.
- 4. When removing the lens stage, DO NOT twist it laterally (Fig. 9).
- 5. Upon removing the lens stage, confirm that there is no abnormality in the IOL. Should a haptic of the IOL become deformed and protrude out of the hole in the back of the injector body as shown in Fig. 10, DO NOT use it.
- 6. After removing the lens stage, DO NOT fill it with additional OVD.
- 7. Push the plunger slowly when it makes contact with the lens.
- 8. DO NOT push the plunger rapidly or intermittently. Push the plunger forward at a smooth constant rate.
- 9. DO NOT pull back the plunger once the plunger has been pushed forward.
- 10. Keep the inlet (*Kowa* mark) upward while pushing the plunger ahead at a constant rate to release the IOL inside the capsular bag.
- 11. If you feel strong resistance or any abnormal feeling during advancement of the plunger, stop using this device.
- 12. During advancement of the plunger, stop using this device if a trailing haptic is out of the groove along the plunger, or if :
 - the leading haptic should twist and extend forward (Fig. 11);
 - the leading haptic should bend and stretch out (Fig. 12);
 - the trailing haptic should extend out of the nozzle (Fig. 13);
 - the plunger should pass above or under the lens optic or bend the optic irregularly (Fig. 14); or
 - the plunger should go too far toward the right or left side (Fig. 15).
- 13. The IOL should be moved forward until the lens optic is rolled and its edges make secure contact. If you stop advancement short of this location, the plunger may get into and be caught in the lens optic and the lens may be damaged when pushing and inserting the IOL afterward.

14. After the IOL is moved into the near tip of the nozzle, move on to the next procedure promptly (within 20 seconds) and place the lens into the eye. Failure of prompt implantation may cause the lens to become stuck inside the nozzle, possibly resulting in damage to the IOL or improper implantation.

15. During the IOL implantation procedure, DO NOT push the lens to the extent of being exposed from the nozzle tip (Fig. 16).

16. The IOL should be released slowly after the nozzle tip has been located just before the centre of the pupillary area.

CALCULATION OF LENS POWER

The surgeon should determine preoperatively the power of lens to be implanted. The A-constant below is presented as a reference value for lens power calculations. When calculating the exact lens power, it is recommended that calculations should be performed individually based on the equipment used and the operating surgeon's experience.

Model	A-constant (A-mode)
YP2.2R / CP2.2R	118.6

Lens power calculation methods are described in the following references.

Sanders, D.R., Retzlaff, J., and Kraff, M.C., "Comparison of SRK II formula and other second generation formulae", *Journal of Cataract and Refractive Surgery*, Vol. 14, pp.136-141, 1988.

Retzlaff, J.A., Sanders, D.R., and Kraff, M.C., "Development of the SRK/T intraocular lens implant power calculation formula", *Journal of Cataract and Refractive Surgery*, Vol. 16, pp.333-340, 1990; *ERRATA*, Vol. 16, pp.528, 1990.

PATIENT IDENTIFICATION CARD

The Patient Identification Card included in the outer packaging should be completed and given to the patient. The patient must provide this card when seeking any future eye care.

HOW SUPPLIED

The Avansee™ Preload1P / Avansee™ Preload1P Clear is supplied sterile in a primary (blister) package for single use only. The primary (blister) package is sterilised with ethylene oxide gas and the package should be opened only under sterile conditions. DO NOT re-sterilise.

EXPIRATION DATE (sterility expiration date)

The sterility expiration date is clearly indicated on the outer package. The date is expressed as four digits for the year and two digits for the month.

RETURN/EXCHANGE POLICY

Regarding lens return or exchange, please contact with your local distributor or Authorized Representative (See below).

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Symbol Explanation:

Symbols	Description	Symbols	Description	Symbols	Description
	Do not re-use / For single use only		Caution		Manufacturer
	Do not re-sterilise		Keep away from sunlight		Authorized representative in the European Community
	Do not use if package is damaged		Keep dry		Dioptre
	Sterilised using ethylene oxide		Temperature limitation		Overall diameter
	Non-sterile		Serial number		Optic diameter (Body diameter)
	Consult instructions for use		Use by (YYYY-MM : year-month)		

Manufacturer:
Kowa Company, Ltd.
 4-14, 3-Chome, Nihonbashi-Honcho, Chuo-ku, Tokyo 103-8433, JAPAN
 Authorized Representative in the European Community:
Kowa Pharmaceutical Europe GmbH
 Görreshof 151, 53347 Alfter, Germany
 Phone: 00-800-2826-7333
 Fax: +49 2222 95891 99
 www.kowapharmaceuticals.eu
 www.avansee.eu

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