





## X Do Not


À ne pas faire / Verbote / Divieti / Nooit / Przeciwwskazania / NA MHN / No




- 9** DO NOT twist laterally when removing the lens stage.  
NE PAS provoquer de torsion latérale lors du retrait du support de stockage de la lentille.  
Linsenhalter beim Entfernen NICHT seitlich drehen.  
NON ruotare lateralmente durante la rimozione della sicura della lente.  
NIET zijwaarts draaien bij verwijdering van de lenspatroon.  
NIE odginać na bok podczas zdejmowania osłony soczewki.  
NA MHN συστρέφετε πλευρικά τη θήκη του φακού κατά την αφαίρεση της.  
NO torcer lateralmente al retirar la plataforma de la lente.




- 10** DO NOT use – should haptic become deformed or protrude from the device.  
NE PAS utiliser – si l'haptique est déformé ou s'il dépasse de l'appareil.  
NICHT verwenden – Haptik verformt oder ragt aus dem Produkt heraus.  
NON utilizzare in caso di deformazione o di protrusione di un'apica dal dispositivo.  
NIET gebruiken – als het haptische deel vervormd is of uit het instrument steekt.  
NIE używać jeżeli część haptyczna uległa odkształceniu lub wystaje z urządzenia.  
MHN χρησιμοποιείτε – όταν το απτικό στοιχείο έχει παραμορφωθεί ή προεξέχει από τη συσκευή.  
NO utilizar – si el háptico se deforma o sobresale del dispositivo.




- 11** DO NOT use – leading haptic twisted and/or extended forward.  
NE PAS utiliser – si l'haptique antérieur est tordu ou s'étend vers l'avant.  
NICHT verwenden – vordere Haptik verdreht und/oder nach vorne gestreckt.  
NON utilizzare in caso di apica di testa ruotata e/o allungata in avanti.  
NIET gebruiken – als het voorste haptische deel gedraaid is en/of naar voren steekt.  
NIE używać – przednia część haptyczna skręcona i wyciągnięta ku przodowi.  
MHN χρησιμοποιείτε – όταν το άνω άκρο του απτικού στοιχείου έχει στρεβλωθεί και/ή εκτείνεται προς τα εμπρός.  
NO utilizar – háptico delantero retorcido y/o extendido.




- 12** DO NOT use – leading haptic bent or stretched out.  
NE PAS utiliser – si l'haptique antérieur est courbé ou dépasse.  
NICHT verwenden – vordere Haptik verbogen oder nach außen gestreckt.  
NON utilizzare – apica di testa piegata o distesa.  
NIET gebruiken – voorste haptisch deel verbogen of gerekt.  
NIE używać: przednia część haptyczna zgięta lub rozciągnięta.  
NA MHN χρησιμοποιείτε – όταν η πρόσθια ακύλη έχει λυγίσει ή εκτείνεται εκτός.  
NO utilizar – háptico delantero doblado o extendido.




- 13** DO NOT use – trailing haptic extended out.  
NE PAS utiliser – si l'haptique postérieur est détendu.  
NICHT verwenden – hintere Haptik herausgestreckt.  
NON utilizzare – apica di coda allungata verso l'esterno.  
NIET gebruiken – achterste haptisch deel steekt uit.  
NIE używać: tylna część haptyczna rozciągnięta.  
NA MHN χρησιμοποιείτε – όταν η οπίσθια ακύλη επεκτείνεται προς τα έξω.  
NO utilizar – háptico trasero extendido.



- 14** DO NOT use – should the plunger pass above or under the lens optic or bend the optic irregularly.  
NE PAS utiliser – si le piston entre dans le corps de la lentille ou s'il tord anormalement la lentille.  
NICHT verwenden – Kolben über oder unter die Linsenoptik oder auf ungewöhnliche Weise in der Optik gebogen.  
NON utilizzare – se lo stantuffo passa sopra o sotto l'ottica della lente o piega l'ottica in modo irregolare.  
NIET gebruiken – indien de zuiger boven of onder de optische lens loopt of de optische lens onjuist buigt.  
NIE używać, jeżeli tłok przechodzi powyżej lub poniżej części optycznej soczewki lub zagina tę część w nieregularny sposób.  
NA MHN χρησιμοποιείτε – όταν το έμβολο έχει περάσει πάνω ή κάτω από την οπτική διάταξη του φακού ή έχει κάμψει ακανόνιστα το φακό.  
NO utilizar – si el émbolo pasara por encima o por debajo de la óptica de la lente o doblara la óptica de manera irregular.



- 15** DO NOT use – plunger has moved too far towards the left or right side.  
NE PAS utiliser – si le piston est allé trop à droite ou à gauche.  
NICHT verwenden – Kolben zu weit nach links oder rechts verschoben.  
NON utilizzare – stantuffo posizionato troppo a destra o a sinistra.  
NIET gebruiken – zuiger te ver naar links of rechts bewegen.  
NIE używać: tłok przesunął się zbyt daleko w lewo lub w prawo.  
NA MHN χρησιμοποιείτε – όταν το έμβολο έχει προχωρήσει πολύ προς τη αριστερή ή τη δεξιά πλευρά.  
NO utilizar – el émbolo se ha movido con demasiada rapidez hacia el lado izquierdo o derecho.



- 16** DO NOT use – lens exposed at nozzle tip before insertion.  
NE PAS utiliser – si la lentille a été exposée à l'extrémité de l'embout avant insertion.  
NICHT verwenden – Linse vor dem Einsetzen auf Injektorspitze gelangt.  
NON utilizzare – lente sporgente sulla punta dell'ugello prima dell'inserimento.  
NIET gebruiken – lens bij punt Injector opening blootgesteld voor het inbrengen.  
NIE używać: soczewka wystaje z końcówki dyszy przed wprowadzeniem.  
NA MHN χρησιμοποιείτε – όταν ο φακός έχει βγει από το άκρο του ρύγχους πριν την εισαγωγή στο μάτι.  
NO utilizar – lente expuesta en la punta de la boquilla antes de la inserción.

## PRESCRIPTION DEVICE

Cauton: This device is restricted to sale by or on the order of an ophthalmologist.

## DEVICE DESCRIPTION

The Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear from Kowa Company, Ltd. is a foldable posterior chamber aspheric intraocular lens (IOL) (Fig.2) that is preloaded in a single-use Kowa original injector (Fig.1). The optic of the lens and the modified-C haptics are made from a UV-absorbing hydrophobic soft acrylic material which, in the case of the Yellow type also contains proprietary blue-light filtering. The spectral transmittance of the IOL (Yellow type) closely replicates that of the natural crystalline lens (Fig.3). Prior to insertion, the optical portion of the lens is folded to allow for placement through a small incision. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance. The anterior surface of the lens is designed as aspheric and the cylinder power is added to posterior surface. Cylinder axis marks indicating flattest meridian are printed on the toric surface.

## INTENDED USE

Kowa's Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This allows the lens to function as a refractive medium in the visual correction of aphakia with corneal astigmatism.

## INDICATIONS

Kowa's Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear IOL is placed in a capsular bag and is designed for implantation after extracapsular cataract extraction or phacoemulsification of cataracts.

## WARNINGS

- Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting the Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear lens in a child or a "special care patient" who has one or more of the following conditions. Before surgery Corneal endothelial damage / Glaucoma / Uveitis / Diabetic retinopathy / Retinal detachment / Congenital ocular anomalies / Choroidal haemorrhage / Shallow anterior chamber / Microphthalmos / Corneal dystrophy / Optic nerve atrophy / Ocular hypertension / Amydriasis / Amblyopia / History of keratoplasty / Iritis / Corneal anomalies / Macular degeneration / Retinal degeneration / Atopic disease / Pseudoexfoliation syndrome and zonular weakness / Zonular rupture and lens luxation (including lens subluxation) / Rubeosis iridis / Patients with intraoperative occurrences of any severe adverse event / or patients who, assessed by a surgeon, require special care for reasons such as an accompanying systemic or ophthalmic disease.
- Prior to surgery, the risks and benefits associated with the implantation of the IOL should be clearly explained by the surgeon to the patient.
- Additional attention is required for "special care patients", including a post-operative follow up by an experienced ophthalmologist with adequate equipment, due to the higher risk of complications and/or insufficient vision recovery. Implantation of the IOL in children should be performed by an ophthalmologist who has sufficient knowledge and experience in paediatric care. This is particularly important when treating a child (under 2 years old) whose small eyeballs can make implantation and the handling of instruments difficult. The chances of IOL exchange are also higher due to the change of eye axial length with the child's growth. Therefore it is important to provide the parents of the young patient with informed consent information prior in any surgery and the IOL should be implanted to the young patients carefully after considering the benefit/risk ratio.
- For active uveitis and a child with uveitis, inflammation should be suppressed by medical treatment prior to intraocular lens implantation, as surgical invasion may cause aggravation of the uveitis or another complication.
- Rotation of the lens away from its intended axis may reduce the astigmatic correction. Misalignment of 30 degrees or greater may deteriorate postoperative astigmatism. If reoperation for axis correction is required, it should be performed as early as possible before the lens adheres to the capsule bag.
- As with any surgical procedure, potential complications accompanying IOL implantation can occur. Adverse events and malfunctions may require discontinuation of the implantation, since they can lead to blindness, permanent problems with visual acuity and grave health hazard or may require IOL extraction and IOL exchange.
- Complications accompanying IOL implantation may include but are not limited to:  
<Adverse events>  
Corneal oedema / Keratitis (including corneal erosion) / Corneal endothelial damage / Acute corneal decompensation / Detachment of

Descemet's membrane / Conjunctivitis and subconjunctival haemorrhage / Hyphaema / Hypopyon / Iris damage / Iritis (iridocyclitis) / Iris adhesion / Iris prolapse / Abnormal pupil (block, capture, deformation, dilatation, etc.) / Uveitis / Zonular rupture / Posterior capsular rupture / After cataract / Vitreous haemorrhage and opacity / Vitreous prolapse / Detachment, hole, tear, etc. of retinal tissue (including macula) / Retinal detachment / Choroidal detachment / Choroid haemorrhage / Macular oedema and degeneration / Expulsive haemorrhage / Endophthalmitis / Fibrin precipitation / Secondary glaucoma / Intraocular pressure elevation (including transient elevated intraocular pressure and ocular hypertension) / Intraocular pressure lowering / Dysphotopsia / Deterioration of visual function (visual acuity and contrast sensitivity) / Error of predicted refractive power / Wound leak.

## <Malfunctions>

Optic damage (breakage, scratch, etc.) / Haptic damage (breakage, scratch, detachment, etc.) / Adhesion of foreign bodies on the lens surface / Lens surface reflection / Lens discoloration or pseudocolouration / Lens opacification (including glistering) / Lens luxation / Lens decentration / Lens dislocation / Lens jamming.

- Back-ups for the intraocular lens and the insertion device should be made available prior to surgery as a part of emergency preparations.
- DO NOT re-use the lens. For single use only.
- DO NOT re-sterilise. Re-use or re-sterilisation of any component may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Re-use or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to the transmission of infectious disease(s) from one patient to another.
- DO NOT store the lens in direct sunlight or in any hot and humid places. Store the lens at 25 °C or lower and keep it dry.

## PRECAUTIONS

- Before using this device, please ensure that you have carefully read and understood this document in order to complete the procedure safely.
- DO NOT use this lens for the patients with capsular rupture or with zonular rupture, or with large posterior capsulotomy.
- Appropriate IOL model and placement of axis should be determined. The location and size of the incision may affect the amount and axis of residual corneal astigmatism. Kowa provides a calculation tool ([www.avansee.eu/](http://www.avansee.eu/)) for surgeons and recommends its use.
- This device should be placed in the posterior chamber. Implantation of the lens in the anterior chamber has not been evaluated for safety and efficacy.
- The lens should not be placed in the ciliary sulcus.
- Before opening the Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear packaging, confirm the description including the IOL model, dioptric and cylinder lens power, and expiration date located on the label.
- The temperatures of the device and the ophthalmic viscosurgical device (OVD) during use should be kept between 21 °C and 25 °C.
- The inside of the aluminium packaging is not sterile. In a non-sterile environment, remove its contents including the primary (blister) package from the aluminium packaging.
- Should any damage or abnormality be found on the primary package after opening the aluminium packaging, DO NOT use, as the device may no longer be sterile.
- After opening the primary package, the device should be handled aseptically.
- When removing the device from the primary package, ensure to be careful not to drop it.
- Before using the device, carefully examine all parts for damage, or any other abnormality. Should any abnormality be found, DO NOT use it.
- If the lens stage is dislodged before use, DO NOT use it.
- DO NOT remove the lens stage before the injection of the OVD.
- DO NOT open the cap on the injector body.
- If there are any abnormalities during the procedure, immediately stop using this device.
- Except for aborting lens insertion, do not interrupt insertion after the lens starts coming out from the nozzle, and keep pushing the plunger until trailing haptic is completely extruded from the nozzle.
- When implantation has been stopped during the procedure due to an abnormality, DO NOT re-use the device and discard it.

## INSTRUCTIONS FOR USE

<Preparation for lens implantation>

- Prior to implanting, examine the lens package for type, power and proper configuration.

















**INSTRUCCIONES DE USO**

<Preparación de la implantación de la lente>

1. Antes de realizar el implante, examine el envase de la lente con respecto al tipo, potencia y configuración adecuados.
2. Abra el envase principal (blíster) y retire el dispositivo en un entorno estéril.
3. Inserte la aguja de inyección para el OVD en la toma de entrada ubicada en la tapa sobre el cuerpo del inyector y llene la boquilla cerca de la línea de puntos mostrada por el diagrama (Fig. 4) con el OVD para inyectar al menos 0,17 ml aproximadamente.

<Implantación de la lente>

4. Sujutando firmemente con la mano la parte de la marca *Kowa* del cuerpo del inyector, sujete la plataforma de la lente entre su pulgar y dedo corazón en ambos lados y sostenga su parte delantera con su dedo índice (Fig. 5), y retire la plataforma de la lente lentamente del cuerpo del inyector recta (Fig. 6).
5. Presione el émbolo a una velocidad lenta y constante; mueva la LIO hacia delante y deténgala en el punto en el que la óptica de la LIO se enrolle y sus bordes entren en contacto de forma segura (Fig. 7). Inmediatamente (en un plazo de 20 segundos) desde este estado, coloque la LIO dentro del ojo.
6. Inserte la punta de la boquilla con la toma de entrada (marca *Kowa*) hacia arriba a través de una incisión justo delante de zona pupilar central (Fig. 8).
7. Mientras empuja el émbolo hacia delante a una velocidad constante, libere la LIO dentro del saco capsular.
8. Presione hasta que el háptico trasero salga fuera completamente y se libere la LIO, saque el cuerpo del inyector de la cámara anterior. A continuación, establezca la LIO dentro del saco capsular utilizando un gancho u otras herramientas.
9. Debe prestarse atención al retirar el OVD del ojo al final de la valoración.

Precauciones en el manejo:

1. Debe inyectarse una cantidad suficiente de OVD para llenar la boquilla y cubrir toda la óptica de la lente insertando la aguja de inyección en la toma de entrada ubicada en la tapa del cuerpo del inyector (hasta el nivel de al menos 0,17 ml desde la toma de entrada). Un OVD insuficiente podría dañar la LIO debido a una lubricación inadecuada dentro de la boquilla.
2. NO deje el inyector con el OVD lleno dentro de la boquilla. Comience inmediatamente el siguiente procedimiento.
3. NO avance el émbolo antes de retirar la plataforma de la lente.
4. Al retirar la plataforma de la lente, NO la tuerza lateralmente (Fig. 9).
5. Tras retirar la plataforma de la lente, confirme que no existe ninguna anomalía en la LIO. Si un háptico de la LIO se deformara y saliera fuera del agujero en la parte posterior del cuerpo del inyector tal y como se muestra en la Fig. 10, NO la utilice.
6. Tras retirar la plataforma de la lente, NO añada OVD adicional.
7. Presione el émbolo lentamente cuando hace contacto con la lente.
8. Empuje el émbolo hacia adelante a una velocidad constante.
9. NO tire hacia atrás del émbolo una vez que lo ha empujado hacia delante.
10. Mantenga la toma de entrada (marca *Kowa*) hacia arriba mientras empuja el émbolo hacia delante a una velocidad constante para liberar la LIO dentro de la bolsa capsular.
11. Si sintiera una fuerte resistencia o una sensación anómala durante el avance del émbolo, deje de utilizar este dispositivo.
12. Durante el avance del émbolo, deje de utilizar este dispositivo si un háptico trasero estuviera fuera de la ranura a lo largo del émbolo, o en caso de que:
  - el háptico delantero se torciera y se extendiera hacia delante (Fig. 11);
  - el háptico delantero se doblase y extendiera (Fig. 12);

- el háptico trasero debería extenderse fuera de la boquilla (Fig. 13);
  - el émbolo pasara por encima o por debajo de la óptica de la lente o doble la óptica de manera irregular (Fig. 14); o
  - el émbolo fuera demasiado rápido hacia la parte derecha o izquierda (Fig. 15).
13. La LIO debe moverse hacia delante hasta que la óptica de la lente se enrolle y sus bordes hagan un contacto seguro. Si detiene el avance en esta ubicación, el émbolo podría entrar dentro y quedarse atrapado en la óptica de la lente y ésta podría resultar dañada al empujar e insertar la LIO posteriormente.
  14. Tras mover la LIO dentro de la punta de la boquilla, pase al siguiente procedimiento inmediatamente (en el plazo de 20 segundos) y coloque la lente en el ojo. Si no se realizara una implantación inmediata, la lente podría quedar adherida dentro de la boquilla, lo que posiblemente podría causar daños en la LIO o una implantación inadecuada.
  15. Durante el procedimiento de implantación de la LIO, NO empuje la lente en la medida en que pueda quedar expuesta fuera de la punta de la boquilla (Fig. 16).
  16. La LIO debe liberarse lentamente tras situar la punta de la boquilla justo delante del centro de la zona pupilar.
  17. Todos los OVD residuales sobre y bajo la lente deben ser retirados cuidadosamente, de lo contrario se podría causar la rotación y alineación incorrecta de la lente.

**CÁLCULO DE LA POTENCIA DE LALENTE**

El cirujano debe determinar antes de la intervención quirúrgica la potencia cilíndrica de la lente que se implantará. La constante A de la siguiente tabla representa un punto de referencia para el cálculo de la potencia de la lente. La constante A y las potencias cilíndricas en el plano LIO inferiores son presentados como valor de referencia para el cálculo de la potencia de la lente. Al calcular la potencia exacta de la lente, se recomienda realizar los cálculos de forma individual en función de los equipos utilizados y de la experiencia del cirujano.

Modelo	Constante A (modo A)
YP-Tx / CP-Tx	118,6

**Potencias cilíndricas (D):**

Modelo: YP-/CP-	T3	T4	T5	T6	T7
Plano LIO (etiquetado)	1,50	2,25	3,00	3,75	4,50
Plano corneal	1,0	1,5	2,0	2,5	3,0

Los métodos de cálculo de la potencia de la lente se describen en las siguientes referencias.

Sanders, D.R., Retzlaff, J., and Kraff, M.C., "Comparison of SRK II formula and other second generation formulae", *Journal of Cataract and Refractive Surgery*, Vol. 14, pp.136-141, 1988.

Retzlaff, J.A., Sanders, D.R., and Kraff, M.C., "Development of the SRK/T intraocular lens implant power calculation formula", *Journal of Cataract and Refractive Surgery*, Vol. 16, pp.333-340, 1990; ERRATA, Vol. 16, pp.528, 1990.

**TARJETA IDENTIFICADORA DE PACIENTES**

La Tarjeta identificadora de pacientes incluida en el envase exterior debe completarse y entregarse al paciente. El paciente debe proporcionar esta tarjeta cuando busque alguna futura asistencia ocular.

**MODO DE SUMINISTRO**

Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear se suministra estéril en un envase (blíster) principal para un solo uso. El envase (blíster) primario se esteriliza con gas de óxido de etileno y el envase debe abrirse solo en condiciones de esterilidad. NO reesterilizar.

**FECHA DE CADUCIDAD (fecha de caducidad de esterilidad)**

La fecha de caducidad de la esterilización se indica claramente en el envase exterior. La fecha se expresa como cuatro dígitos para el año y dos dígitos para el mes.

**POLÍTICA DE DEVOLUCIÓN/CAMBIO DEL ARTÍCULO**

En relación con la devolución o cambio de la lente, póngase en contacto con su distribuidor local o representante autorizado (véase abajo).

# Avansee™ Preload1P Toric

# Avansee™ Preload1P Toric Clear

YP-Tx / CP-Tx (x: 3, 4, 5, 6 or 7)

Rev. 001  
2019 - 01

Package insert – Instructions for use

Feuille de l'emballage – Manuel d'utilisation

Packungsbeilage – Gebrauchsanleitung

Foglietto illustrativo – Istruzioni per l'uso

Bijlsuiter – Gebruiksaanwijzing

Ulotka dołączana do opakowania – Instrukcja użycia

Ένθετο στη συσκευασία – Οδηγίες χρήσης

Prospecto – Instrucciones de uso

Explicación de los símbolos:

Símbolos	Descripción	Símbolos	Descripción	Símbolos	Descripción
	No reutilizar / para un solo uso		Precaución		Fabricante
	No reesterilizar		Mantener fuera de la luz solar		Representante autorizado para la Comunidad Europea
	No utilizar si el envase está dañado		Mantener en lugar seco	D	Dioptrías
	Esterilizado con óxido de etileno		Limitación de temperatura	φ T	Diámetro total
	No estéril		Número de serie	φ B	Diámetro óptico (diámetro del cuerpo)
	Consultar las instrucciones de uso		Fecha de caducidad (AAAA-MM: año-mes)		

Fabricante:

*Kowa Company, Ltd.*  
4-14, 3-Chome, Nihonbashi-Honcho, Chuo-ku, Tokyo 103-8433, JAPÓN

Representante autorizado para la Comunidad Europea:

*Kowa Pharmaceutical Europe GmbH*  
Görreshof 151, 53347 Alfter, Alemania  
Tel.: 00-800-2826-7333  
Fax: +49 2222 95891 99  
www.kowapharmaceuticals.eu      www.avansee.eu



CE0197

